

# Eleganza Sito



Transport and rest chair

# Eleganza Sito – comfort and functionality

The Eleganza Sito transport and rest chairs bring comfort to the patient and, thanks to easy manipulation and the variety of adjustments, also help the staff.

- ERGONOMIC DESIGN
- BROAD RANGE OF ADJUSTMENT
- EXCELLENT MANOEUVRABILITY
- VARIOUS ACCESSORIES



model D

▪ Removable head rest

▪ Adjustable back rest

▪ Adjustable and removable arm rest (6 positions)

▪ Adjustable seat (seat assistance)

▪ Central braking system with lockable castor for steering

▪ Adjustable foot rest (can be hidden)

▪ Height adjustment



model A



model B

## TECHNICAL PARAMETERS

model	R1	R2
Width	61 cm	66 cm
Length	74 cm	74 cm
Height	112 cm	112 cm
Weight	21 kg	22 kg
Max. user weight	135 kg	

model	A	B
Width	67 cm	68 cm
Length	79 cm	76 cm
Height	112 cm	119 cm
Weight	28 kg	32 kg
Max. user weight	135 kg	135 kg

model	C/D
Width	76 cm
Length	95 cm
Min. height	125 cm
Max. height	154 cm
Back rest adjustment	12°–45°
Weight	58 kg
Max. user weight	135 kg

### BENEFIT FOR THE PATIENT

- The chair can be used for partial verticalisation which increases the mobilisation of the patient.
- User-friendly, comfortable chair.
- Unique ergonomic design of sitting area and back rest.
- Broad range of adjustments allows position for maximum comfort to be achieved.

### BENEFIT FOR THE STAFF

- Multi-purpose solution allows personnel to choose optimal model.
- Various functions make the daily routine easier and safer.
- Range of accessories optimizes the chair according to the health condition of the patient (IV poles, pair of side supports, small body kit, eating board, leg support...).