

Eleganza **Mano**



Patient bedside cabinet

Eleganza Mano – easy manipulation and operation

The unique construction of the Eleganza Mano patient bedside cabinet has been designed with patient independence, infection control and ease of handling in mind.

- PATIENT INDEPENDENCE
- SAFE ROUNDED EDGES
- EASY MANIPULATION
- EASY CLEANING



EASE OF HANDLING

- Combination of central braking system, handrails and position of castors ensure that the patient can easily manipulate the bedside cabinet.
- Ergonomically designed top and bottom drawers can be accessed from both sides allowing ease of manipulation and patient access.
- Patient can independently access eating board, towel holder, shelf for bedpan and other storage spaces.

INFECTION CONTROL

- Easily accessible central braking mechanism can be activated by the patient using just one hand, allowing the bedside cabinet to be moved to a convenient position.
- Smooth and non-absorbent surfaces ensure that the bedside cabinet is easy to clean and decontaminate.
- All drawers and internal trays can be removed to enable thorough cleaning of the cabinet.



Easily accessible central braking control



Drawers can be accessed from both sides allowing easy patient access.

TECHNICAL DATA

| | |
|-------------------------------|--------------|
| External dimensions: | 51 x 48,5 cm |
| Upper desk dimensions | 380x 41 cm |
| Height (incl. gallery) | 93 cm |
| Height of upper desk | 89 cm |
| Eating board dimensions | 60 x 35 cm |
| Height of dining board | 73 – 110 cm |
| Load capacity of eating board | 5 kg |
| Castors | 7,5 cm |